



# A robust fetal ECG detection method for abdominal recordings

Porcelli Angelica – 78083A

Biomedical Signal Processing



# Table of content

- ▶ Introduction
- ▶ The algorithm
- ▶ Dataset
- ▶ Result
- ▶ Conclusion



# Fetal ECG Extraction Challenges

When recording **abdominal ECG** from a *pregnant woman*, the signal is a **mixture** of different physiological and non-physiological components:


## 1. Maternal ECG (MECG)



## 2. Fetal ECG (FECG)



## 3. Noise and artifacts







# What is the paper approach ?

The paper proposes a **sequential, physiologically-informed approach**.

The idea is:

instead of trying to separate everything at once (as ICA does), remove problems **one by one**, in a meaningful order.

The algorithm:

**7 step**, and **each block removes a specific distortion**.






# Table of content

- ▶ Introduction
- ▶ The algorithm
- ▶ Dataset
- ▶ Result
- ▶ Conclusion



# The Algorithm

1. Remove baseline
2. Remove powerline
3. Increase samplings
4. Maternal QRS detector
5. MECG Canceller
6. Fetal QRS detector
7. FECG detector








# 1. Remove Baseline

Objective: **remove slow oscillations** caused by breathing, body movements, ...

High-pass FIR filter :

- 1000 + 1 coefficients
- 3 Hz cutoff frequency
- Window Hamming





$S_1 \rightarrow S_2$



## 2. Remove Powerline

Objective: **remove 50 Hz electrical interference** and its **harmonics**



Use of **adaptive noise cancelling (ANC)**:

- Estimates the 50 Hz interference using *model*

$$\hat{d}(t) = A(t)\sin(\phi(t)) + B(t)\cos(\phi(t))$$


Clean Signal → Update of Coefficients → Phase Adaptation  
using an LMS rule driven  
by the residual error

$S_2 \rightarrow S_3$





# 3. Increase samplings

Objective: **improve temporal resolution** and QRS alignment



The signal is upsampled to 2000 Hz.

- More accurate MECG canceller.
- More visible FECG.



```
resample_poly(X, L, M, axis=0, window=('kaiser', 5.0))
```

$S_3 \rightarrow S_4$



# 4. M-QRS Detection

Objective: **identify maternal QRS** complexes to enable correct ECG cancellation

## 1. PCA enhancement

PCA-based multichannel QRS enhancement

Extract PC1

## 2. Rough detection




Cross-correlation with a synthetic QRS template

## 3. Real Template

Averaging of  $\pm 50$  ms segments around subset

## 4. Final detection



Cross-correlation with QRS template



$S_4 \rightarrow \text{M-QRS}$



# 4. M-QRS Detection




refractory\_ms = 300, refine\_ms =  $\pm 40$ , duration\_ms = 100

$S_4 \rightarrow \text{M-QRS}$



# 5. MECG Canceller

Objective: remove maternal ECG (MECG) while preserving fetal ECG



## 1. Extract window

$$m(t) = S4[R - 250\text{ ms} : R + 450\text{ ms}]$$

## 2. Create a template

$$\mu(t) = \frac{1}{N} \sum_{i=1}^N m_i(t)$$

## 4. Find adaptive scaling


$$\underline{a} = (M^T M)^{-1} M^T \underline{m}$$

$$\hat{\underline{m}} = M \underline{a}$$

## 5. Template subtraction

$$S5(t) = S4(t) - \hat{m}(t)$$


## 3. Split and create matrix



S4 → S5




# 5. MECG Canceller




```
for i in range(3):
    S4 = mecg_canceller(S4, mQRS, fs_up, N=N)
```

Zoom around maternal QRS after cancellation



AECG4



S4 → S5





# 6. F-QRS Detection

Objective: identify fetal QRS complexes

1. PCA enhancement
2. Rough detection
3. Real Template
4. Final detection



short refractory period  
short peak refinement  
short duration




$S_5 \rightarrow \text{fQRS}$



# 6. F-QRS Detection



`refractory_ms = 150,`  
`refine_ms = ±20,`  
`duration_ms = 40`




S5 → fQRS



# 7.1 FECG Extraction

Objective: obtain a clean and reliable estimate of the fetal **cardiac waveform**



## 1. Window extraction:


For each fetal QRS, a  $\pm 100$  ms window centered on the R-peak.

## 2. Intrinsic alignment:

Since each window is centered on the detected QRS, all segments are aligned.

## 3. Ensemble averaging:

The first 150\* aligned segments are averaged, reducing uncorrelated noise by a factor of  $\sqrt{150} \approx 12$ .




S5 → S6





# 7.2 FHR

Objective: **evaluate the temporal consistency** and correctness of fetal QRS detection

$$RR_i = \frac{R_{i+1} - R_i}{f_s}$$



$$FHR_i = \frac{60}{RR_i}$$



Mean bpm 130.21 detected  
Mean bpm 128.76 real

S5 → S6





# Table of content

- ▶ Introduction
- ▶ The algorithm
- ▶ Dataset
- ▶ Result
- ▶ Conclusion



# Dataset

PhysioNet Challenge 2013 (Set A): <https://physionet.org/content/challenge-2013/1.0.0/>

## Characteristics:

25 real-world recordings


4 simultaneous abdominal channels

Duration: 60 sec

Sampling frequency: 1000 Hz

*+ fQRS annotations*

## AECG1, AECG2, AECG3, AECG4







# Table of content

- ▶ Introduction
- ▶ The algorithm
- ▶ Dataset
- ▶ **Result**
- ▶ Conclusion




# Signal quality metrics

$$Reliability = 1 - \frac{num\_outliers}{total\_point\_FHR}$$


$$num\_outliers = |FHR_i - median| > 10 \text{ bpm}$$

$$SUCCESS = (60 \leq mean\_bpm \leq 220) \text{ and } (reliability > 0.60)$$

Comparison Across Groups





Reliability Record







# Signal quality metrics





$$SNR_j = \frac{P_{F_j}}{P_{N_j}} = \frac{P_{S_{6,j}}}{P_{S_{5,j}-S_{6,j}}},$$

$$SIR_j = \frac{P_{F_j}}{P_{M_j}} = \frac{P_{S_{6,j}}}{P_{S_{4,j}-S_{5,j}}}.$$










# Comparison with FastICA





# Comparison with FastICA







# Table of content


- ▶ Introduction
- ▶ The algorithm
- ▶ Dataset
- ▶ Result
- ▶ Conclusion



# Conclusion

**Non-blind methods outperform blind ICA.**

The **paper-pipeline** is more **controllable** and consistent with **SNR/SIR**, while **FastICA** shows **higher variability**.







# Thanks

Porcelli Angelica – 78083A  
Biomedical Signal Processing